

# Split-Specimen Cup™

Innovative drug test cup design eliminates urine handling and donor tampering, providing a unique and efficient approach for on-site drug testing.

The Peeler Group  
The Peeler Group

## Test With Confidence.

- Split-Specimen Collection
- Screw Top Lid for Secure Closure
- Collector Activated vs Donor Activated
- Easy to Photocopy Results
- Built-in Adulteration Testing Available
- FDA Approved




For More Information  
888.842.8274


Point of care medical diagnostic testing products

For More Information  
888.842.8274

**1 Collect Specimen**




**2 Insert Key**



**3 Peel Back Label**



**4 Read the Results**



Read results at (5) minutes.  
Results stable up to (60) minutes.

All positive results are presumptive and should be confirmed by an alternate method. (e.g. GC/MS)

Drug Name:	Abbreviation:	Cutoff Levels:
Amphetamine	AMP	1000ng/mL
Barbiturates	BAR	300ng/mL
Benzodiazepines	BZO	300ng/mL
Buprenorphine	BUP	10ng/mL
Cocaine	COC	300ng/mL
Cocaine	COC150	150ng/mL
Ecstasy	MDMA	500ng/mL
Marijuana	THC	50ng/mL
Methadone	MTD	300ng/mL
Methamphetamine	mAMP	1000ng/mL
Morphine	MOR	300ng/mL
Opiates	OPI	2000ng/mL
Oxycodone	OXY	100ng/mL
Phencyclidine	PCP	25ng/mL
Propoxyphene	PPX	300ng/mL
Tricyclic Antidepressants	TCA	1000ng/mL

**GCMS Confirmation Services Available**

